



Additional Stirrer Blades for Bola Stirrer Shafts, PTFE

Additional Stirrer Blades for Bola Stirrer Shafts, PTFE

Useful extension to the existing range of Bola stirrer shafts. The additional blades can optionally be fixed in your preferred position on every stirrer shaft with diameter 8, 10 or 16 mm (a flat spanner is included). The stirrer blades are entirely made of PTFE, the mounting devices are made of a special PTFE-compound. A high chemical and thermal (-200 to +250 °C) resistance is assured. This "building block system" allows the user to define a stirrer of their own special requirements.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG06227770 | Stirrer Blades for Bola Stirrer Shafts, PTFE, Maxi propeller-three-bladed, Shaft diam. 16 mm | |
| LLG06228187 | Solo stirrer shaft 600 mm, Ø shaft 10 mm | |
| LLG06233277 | Additional Stirrer Blades for Bola Stirrer Shafts, PTFE, Description Propeller-three-bladed , Shaft diam. 8 mm, Length 18 mm, Paddle diam. 75 mm | 1pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

