

**pH combination electrode LE 438**

Mettler-Toledo Online

pH combination electrode LE 438

With epoxy body and gel electrolyte, fixed cable with BNC / Cinch connection.
For simple measurements, particularly suitable for portable pH meters, with temperature sensor.

Specifications

Electrolyte: 3 mol / L KCl sat. with AgCl as gel, non-refillable

Body length: 120 mm

pH: 0 ... 14 pH

Temperature: 0 ... 80 °C

Temperature sensor: NTC 30 kW



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG09041400 | pH combination electrode LE 438, Type LE 438, Connection BNC / Cinch | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

