

**Thermohygrometer, Klimatherm****TFA Dostmann****Thermohygrometer, Klimatherm**

Spun metal case with dial, separate humidity and temperature indicators, and brass hanging ring.

**Measuring range**

**Relative humidity: 15 to 100 %**

**Temperature: 0 to 40 °C**

**Weight: 61 g**



| Code        | Description                                | Packaging |
|-------------|--------------------------------------------|-----------|
| LLG09726404 | Thermohygrometer, Klimatherm, Diam. 102 mm | 1 pz      |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

