



Temperature sensor for RET control-visc

IKA

Temperature sensor for RET® control-visc



Code	Description	Packaging
LLG06270208	Temperature sensor for RET® control-visc, Type PT 1000.70, Stainless steel,, Immersion depth 230 mm, Diam. 3 mm	
LLG09720348	Temperature sensor for RET® control-visc, Type PT 1000.50, Stainless steel, dual probe,, Immersion depth 230 mm, Diam. 3 mm	1 pz.
LLG09720349	Temperature sensor for RET® control-visc, Type PT 1000.51, Stainless steel, glass coated, *, Immersion depth 230 mm, Diam. 8 mm	1 pz.
LLG09720351	Temperature sensor for RET® control-visc, Type PT 100.51, Stainless steel, glass coated, *, Immersion depth 230 mm, Diam. 8 mm	1 pz.
LLG09720352	Temperature probe PT 1000.80	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

