



## Overgoggles uvex super OTG 9169

### Overgoggles uvex super OTG 9169

The super OTG and the super f OTG are equipped with this state-of-the-art lens technology.

Alongside the tried and tested hingeless spectacles and highly flexible soft ear pieces, a hinged design has been developed.

The revolutionary new design of the super f OTG with hinge combines all the principles of contemporary spectacles with the advantages of state-of-the-art materials.

Both designs are well-suited to any wearer situation, and can also be used as overspecs for prescription wearers.

#### Coating of the lens

1. uvex supravision sapphire: excellent scratch and chemical resistance, 100% UV protection
2. uvex supravision excellence: excellent scratch and chemical resistance and easy to clean ( nano clean hybrid )



Code	Description	Packaging
LLG07653307	Protection glasses super OTG 9169colour	1 pz.
LLG06234320	Overgoggles uvex super OTG 9169, Type super OTG, Colour navy blue, Lens PC transparent/UV 2-1.2, Coating 2	1 pz.
LLG07650112	Overgoggles uvex super OTG 9169, Type super OTG, Colour navy blue, Lens PC transparent/UV 2-1.2, Coating 1	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

