



Overgoggles uvex super OTG 9169

Overgoggles uvex super OTG 9169

The super OTG and the super f OTG are equipped with this state-of-the-art lens technology.

Alongside the tried and tested hingeless spectacles and highly flexible soft ear pieces, a hinged design has been developed.

The revolutionary new design of the super f OTG with hinge combines all the principles of contemporary spectacles with the advantages of state-of-the-art materials.

Both designs are well-suited to any wearer situation, and can also be used as overspecs for prescription wearers.

Coating of the lens

1. uvex supravision sapphire: excellent scratch and chemical resistance, 100% UV protection
2. uvex supravision excellence: excellent scratch and chemical resistance and easy to clean (nano clean hybrid)



Code	Description	Packaging
LLG07653307	Protection glasses super OTG 9169colour	1 pz.
LLG06234320	Overgoggles uvex super OTG 9169, Type super OTG, Colour navy blue, Lens PC transparent/UV 2-1.2, Coating 2	1 pz.
LLG07650112	Overgoggles uvex super OTG 9169, Type super OTG, Colour navy blue, Lens PC transparent/UV 2-1.2, Coating 1	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

