



Swiss Boy laboratory jacks, anodised aluminium.

Swiss Boy laboratory jacks, anodised aluminium.

Anodised aluminium. Corrosion proof, stable and easily adjustable. Upper and lower end positions are adjustable. At reaching the setting, the laboratory support stops automatically. In various colours as outlined below.



Code	Description	Packaging
LLG09117110	Swiss Boy laboratory jacks, anodised aluminium., Colour Pink, Width 140 mm, Depth 120 mm, Min. height 60 mm, Max. height 243 mm, Max. stat. 30 kg	1 pz.
LLG09117115	Swiss Boy laboratory jacks, anodised aluminium., Colour Gold, Width 150 mm, Depth 122 mm, Min. height 80 mm, Max. height 278 mm, Max. stat. 40 kg	1 pz.
LLG09117100	Swiss Boy laboratory jacks, anodised aluminium., Colour Green, Width 40 mm, Depth 50 mm, Min. height 31 mm, Max. height 93 mm, Max. stat. 10 kg	1 pz.
LLG09117118	Swiss Boy laboratory jacks, anodised aluminium., Colour Grey, Width 320 mm, Depth 240 mm, Min. height 120 mm, Max. height 500 mm, Max. stat. 70 kg	1 pz.
LLG09117116	Swiss Boy laboratory jacks, anodised aluminium., Colour Light blue, Width 216 mm, Depth 180 mm, Min. height 95 mm, Max. height 398 mm, Max. stat. 60 kg	1 pz.
LLG09117105	Swiss Boy laboratory jacks, anodised aluminium., Colour Blue, Width 75 mm, Depth 80 mm, Min. height 49 mm, Max. height 147 mm, Max. stat. 15 kg	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

