




120 ° Y connectors, PP

Kartell

120 ° Y connectors, PP
Autoclavable.



Code	Description	Packaging
LLG09207231	120 ° Y connectors, PP, Diam. mm 4/5/6 min./max., Tubing i.d. 1.5 mm	1 pz
LLG09207232	120 ° Y connectors, PP, Diam. mm 6/7/8 min./max., Tubing i.d. 2.7 mm	1 pz
LLG09207233	120 ° Y connectors, PP, Diam. mm 8/9/10 min./max., Tubing i.d. 4.7 mm	1 pz
LLG09207234	120 ° Y connectors, PP, Diam. mm 10/11/12 min./max., Tubing i.d. 5.5 mm	1 pz
LLG09207235	120 ° Y connectors, PP, Diam. mm 12/13/14 min./max., Tubing i.d. 7.6 mm	1 pz
LLG09207236	120 ° Y connectors, PP, Diam. mm 14/15/16 min./max., Tubing i.d. 9.6 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

