

**Magnetic stirrer, operated by water/air pressure****Kartell****Magnetic stirrer, operated by water/air pressure**

PP-ABS construction. This unique stirrer is operated either by water pressure between 1kg/cm² and 2.5kg/cm² (15psi to 40psi) or by air pressure. No electrical supply is needed, thus eliminating heat, electric shock and explosion hazards. Corrosion-proof, all-plastic construction. Can be used on the bench or clamped to a stand with a 12mm (½") diameter rod using the integral bosshead.

Specifications

Dimensions (Ø x H): 129 x 61 mm



Code	Description	Packaging
LLG09303043	Magnetic stirrer PP-ABS with air or water pressure Ø 129 mm, height 61 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

