



Test tubes, Borosilicate glass 3.3

Scherf

Test tubes, Borosilicate glass 3.3



Code	Description	Packaging
LLG04656051	Test tubes, Borosilicate glass 3.3, Diam. 30 mm	50 pz.
LLG04658113	Test tubes, Borosilicate glass 3.3, Diam. 10 mm	100 pz.
LLG04658114	Test tubes, Borosilicate glass 3.3, Diam. 12 mm	100 pz.
LLG04658129	Test tubes, Borosilicate glass 3.3, Diam. 20 mm	100 pz.
LLG04658130	Test tubes, Borosilicate glass 3.3, Diam. 25 mm	50 pz.
LLG04658131	Test tubes, Borosilicate glass 3.3, Diam. 16 mm	100 pz.
LLG04658132	Test tubes, Borosilicate glass 3.3, Diam. 18 mm	100 pz.
LLG04658133	Test tubes, Borosilicate glass 3.3, Diam. 20 mm	100 pz.
LLG04658134	Test tubes, Borosilicate glass 3.3, Diam. 25 mm	50 pz.
LLG04658135	Test tubes, Borosilicate glass 3.3, Diam. 30 mm	50 pz.
LLG06281479	Test tubes, Borosilicate glass 3.3, Diam. 25 mm	50 pz.
LLG06281701	Test tubes, Borosilicate glass 3.3, Diam. 12 mm	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

