



Accessories for blender BagSystem

Accessories for blender BagSystem®



Code	Description	Packaging
LLG09570088	Accessories for blender BagSystem®, Type BagClip® 80/100, bag closure clip	50 pz.
LLG09570089	Accessories for blender BagSystem®, Type BagClip® 400, bag closure clip	50 pz.
LLG09570090	Accessories for blender BagSystem®, Type BagClip® 3500, bag closure clip	50 pz.
LLG09570083	Accessories for blender BagSystem®, Type BagOpen® 80/100, bag opener	1 pz.
LLG09570084	Accessories for blender BagSystem®, Type BagOpen® 400, bag opener	1 pz.
LLG09570085	Accessories for blender BagSystem®, Type BagOpen® 3500, bag opener	1 pz.
LLG09570087	Accessories for blender BagSystem®, Type BagRack® 400, bag storage for 10 bags	1 pz.
LLG09570086	Accessories for blender BagSystem®, Type BagRack® 80/100, bag storage for 11 bags	1 pz.
LLG09570091	Accessories for blender BagSystem®, Type BagPipet®, pipette with 3 fixed graduations: 0.1, 0.9, 1.0ml	1 pz.
LLG09570092	Accessories for blender BagSystem®, Type BagTips® Regular, straws for BagPipet®, 19cm, sterile	1000 pz.
LLG09570093	Accessories for blender BagSystem®, Type BagTips® Jumbo, straws for BagPipet®, 24cm, sterile	1000 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

