



## Glass-Metal Syringes SANITEX

### Glass-Metal Syringes SANITEX®

Made from neutral glass, with steel piston. Brown graduation, acid and alkali resistant. With centered Luer-slip/Luer-lock. May only be used for technical purposes. Autoclavable up to 200 °C.



Code	Description	Packaging
LLG09222222	Glass-metal syringes, cap. 2 ml, Luer con	1 pz.
LLG09222232	Glass-Metal Syringes SANITEX®, 20 ml, Grad. 1 ml, Connection Luer-Slip	1 pz
LLG09222242	Glass-metal syringes, cap. 2 ml, Luer-loc	1 pz.
LLG09222250	Glass-Metal Syringes SANITEX®, 10 ml, Grad. 0.5 ml, Connection Luer-Lock	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

