


**Tube racks, PP****Kartell****Tube racks, PP**

Will not float. Alphanumeric grid referenced on top tier for easy identification of samples. Stackable when empty. Can be used in freezers or water baths. Autoclavable at 121°C for 20 min.



Code	Description	Packaging
LLG06204134	Tube racks, PP, For tubes diam. 30 mm, Colour yellow, Dimensions (l x w x h) 300 x 112 x 85 mm, Array 24	1 pz
LLG06207191	Tube racks, PP, For tubes diam. 30 mm, Colour white, Dimensions (l x w x h) 300 x 112 x 85 mm, Array 24	1 pz
LLG06401582	Tube racks, PP, For tubes diam. 30 mm, Colour red, Dimensions (l x w x h) 300 x 112 x 85 mm, Array 24	1 pz
LLG06205110	Tube racks, PP, For tubes diam. 20 mm, Colour white, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 40	1 pz
LLG06206206	Tube racks, PP, For tubes diam. 20 mm, Colour blue, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 40	1 pz
LLG06224461	Tube racks, PP, For tubes diam. 20 mm, Colour yellow, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 40	1 pz
LLG06236328	Tube racks, PP, For tubes diam. 25 mm, Colour white, Dimensions (l x w x h) 295 x 125 x 85 mm, Array 40	1 pz
LLG06236329	Tube racks, PP, For tubes diam. 25 mm, Colour blue, Dimensions (l x w x h) 295 x 125 x 85 mm, Array 40	1 pz
LLG06236330	Tube racks, PP, For tubes diam. 25 mm, Colour yellow, Dimensions (l x w x h) 295 x 125 x 85 mm, Array 40	1 pz
LLG09193201	Tube racks, PP, For tubes diam. 20 mm, Colour red, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 40	1 pz
LLG09193202	Tube racks, PP, For tubes diam. 25 mm, Colour red, Dimensions (l x w x h) 295 x 125 x 85 mm, Array 40	1 pz
LLG06204936	Tube racks, PP, For tubes diam. 16 mm, Colour red, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 60	1 pz
LLG06205109	Tube racks, PP, For tubes diam. 16 mm, Colour white, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 60	1 pz
LLG06206305	Tube racks, PP, For tubes diam. 16 mm, Colour yellow, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 60	1 pz
LLG06206347	Tube racks, PP, For tubes diam. 16 mm, Colour blue, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 60	1 pz
LLG06205483	Tube racks, PP, For tubes diam. 13 mm, Colour white, Dimensions (l x w x h) 246 x 105 x 64 mm, Array 90	1 pz
LLG06206439	Tube racks, PP, For tubes diam. 13 mm, Colour yellow, Dimensions (l x w x h) 246 x 105 x 64 mm, Array 90	1 pz
LLG06224460	Tube racks, PP, For tubes diam. 13 mm, Colour blue, Dimensions (l x w x h) 246 x 105 x 64 mm, Array 90	1 pz



Code	Description	Packaging
LLG09193200	Tube racks, PP, For tubes diam. 13 mm, Colour red, Dimensions (l x w x h) 246 x 105 x 64 mm, Array 90	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

