


**Tube racks, PP****Kartell****Tube racks, PP**

Will not float. Alphanumeric grid referenced on top tier for easy identification of samples. Stackable when empty. Can be used in freezers or water baths. Autoclavable at 121°C for 20 min.



| Code        | Description                                                                                              | Packaging |
|-------------|----------------------------------------------------------------------------------------------------------|-----------|
| LLG06204134 | Tube racks, PP, For tubes diam. 30 mm, Colour yellow, Dimensions (l x w x h) 300 x 112 x 85 mm, Array 24 | 1 pz      |
| LLG06207191 | Tube racks, PP, For tubes diam. 30 mm, Colour white, Dimensions (l x w x h) 300 x 112 x 85 mm, Array 24  | 1 pz      |
| LLG06401582 | Tube racks, PP, For tubes diam. 30 mm, Colour red, Dimensions (l x w x h) 300 x 112 x 85 mm, Array 24    | 1 pz      |
| LLG06205110 | Tube racks, PP, For tubes diam. 20 mm, Colour white, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 40  | 1 pz      |
| LLG06206206 | Tube racks, PP, For tubes diam. 20 mm, Colour blue, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 40   | 1 pz      |
| LLG06224461 | Tube racks, PP, For tubes diam. 20 mm, Colour yellow, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 40 | 1 pz      |
| LLG06236328 | Tube racks, PP, For tubes diam. 25 mm, Colour white, Dimensions (l x w x h) 295 x 125 x 85 mm, Array 40  | 1 pz      |
| LLG06236329 | Tube racks, PP, For tubes diam. 25 mm, Colour blue, Dimensions (l x w x h) 295 x 125 x 85 mm, Array 40   | 1 pz      |
| LLG06236330 | Tube racks, PP, For tubes diam. 25 mm, Colour yellow, Dimensions (l x w x h) 295 x 125 x 85 mm, Array 40 | 1 pz      |
| LLG09193201 | Tube racks, PP, For tubes diam. 20 mm, Colour red, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 40    | 1 pz      |
| LLG09193202 | Tube racks, PP, For tubes diam. 25 mm, Colour red, Dimensions (l x w x h) 295 x 125 x 85 mm, Array 40    | 1 pz      |
| LLG06204936 | Tube racks, PP, For tubes diam. 16 mm, Colour red, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 60    | 1 pz      |
| LLG06205109 | Tube racks, PP, For tubes diam. 16 mm, Colour white, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 60  | 1 pz      |
| LLG06206305 | Tube racks, PP, For tubes diam. 16 mm, Colour yellow, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 60 | 1 pz      |
| LLG06206347 | Tube racks, PP, For tubes diam. 16 mm, Colour blue, Dimensions (l x w x h) 246 x 105 x 72 mm, Array 60   | 1 pz      |
| LLG06205483 | Tube racks, PP, For tubes diam. 13 mm, Colour white, Dimensions (l x w x h) 246 x 105 x 64 mm, Array 90  | 1 pz      |
| LLG06206439 | Tube racks, PP, For tubes diam. 13 mm, Colour yellow, Dimensions (l x w x h) 246 x 105 x 64 mm, Array 90 | 1 pz      |
| LLG06224460 | Tube racks, PP, For tubes diam. 13 mm, Colour blue, Dimensions (l x w x h) 246 x 105 x 64 mm, Array 90   | 1 pz      |



| Code        | Description                                                                                           | Packaging |
|-------------|-------------------------------------------------------------------------------------------------------|-----------|
| LLG09193200 | Tube racks, PP, For tubes diam. 13 mm, Colour red, Dimensions (l x w x h) 246 x 105 x 64 mm, Array 90 | 1 pz      |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

