



Precision refractometers, basic model

Precision refractometers, basic model

Simple and precise optical instruments to analyse liquids and their solid contents. Models are equipped with an automatic temperature compensation device (ATC). Supplied with a plastic box and comfortable leather-type case.



Code	Description	Packaging
LLG09801242	Precision refractometers, basic model, Type RLC/ATC, For	1 pz.
LLG09801200	Precision refractometers, basic model, Type RCZ/ATC, For	1 pz.
LLG09801201	Precision refractometers, basic model, Type RTS, For	1 pz.
LLG09801250	Precision refractometers, basic model, Type RCS/ATC, For Salinity	1 pz.
LLG09801202	Precision refractometers, basic model, Type RPI/ATC, For Biological fluids: urine, serum	1 pz.
LLG09801253	Precision refractometers, basic model, Type RMB/ATC, For Wine	1 pz.
LLG09801255	Precision refractometers, basic model, Type RBB/ATC, For Wine	1 pz.
LLG09801252	Precision refractometers, basic model, Type RMA/ATC, For Wine	1 pz.
LLG09801245	Precision refractometers, basic model, Type RPD/ATC, For	1 pz.
LLG09801246	Precision refractometers, basic model, Type RPU, For	1 pz.
LLG09801254	Precision refractometers, basic model, Type RBO/ATC, For Wine	1 pz.
LLG09801256	Precision refractometers, basic model, Type RGA/ATC, For Alcohol	1 pz.
LLG09801241	Precision refractometers, basic model, Type RBC/ATC, For	1 pz.
LLG09801249	Precision refractometers, basic model, Type RSM/ATC, For Salinity	1 pz.
LLG09801244	Precision refractometers, basic model, Type RFC/ATC, For	1 pz.
LLG09801247	Precision refractometers, basic model, Type RMM/ATC, For Honey and jam	1 pz.
LLG09801251	Precision refractometers, basic model, Type RPC/ATC, For Freezing point batteries	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

