

**Polarimeter PL1****Exacta Optech****Polarimeter PL1**

Simple and sturdy instrument to measure the rotation angle of optically active substances. The measure in 180° is directly readable on two mechanical dials.

Light source: LED lamp (589.3 nm)

Polarimetric Tubes one 100 mm and one 200 mm

Dust cover

Optical Glass Windows (4 pcs)

Rubber seals

Fuses

Range:

Reading precision: 0 to 180°

0.05°



Code	Description	Packaging
LLG09801259	Polarimeter PL1, Type PL1	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

