

**Metal (ORP) electrodes ScienceLine Pt 6280****SI Analytics****Metal (ORP) electrodes ScienceLine Pt 6280**

With glass shaft, pH glass membrane reference system and plug head.

Specifications

Temperature range: -5 ... 100 °C

Length / diam.: 120 mm / 12 mm

Electrolyte: KCl 3 mol/l

Diaphragm: Ceramic

Sensor: Platinum pin (diam. 1 mm)



Code	Description	Packaging
LLG09041211	Metal (ORP) electrodes ScienceLine Ag 6280, Type Ag 6280	1 pz.
LLG09041207	Metal (ORP) electrodes ScienceLine Pt 6280, Type Pt 6280	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

