



Accessories for Reverse osmosis systems Ultra Clear RO / RO EDI / LaboStar RO DI

Accessories for Reverse osmosis systems Ultra Clearä RO / RO EDI / LaboStarä RO DI



Code	Description	Packaging
LLG06259821	Sterile filter 0.2 µm with endotoxin retention, For LaboStar RO DI, Change interval every 6 months	
LLG09914524	Pre-treatment module AMB, For Ultra Clear RO, LaboStar RO/DI, Change interval every 6 months*	1 pz
LLG06285641	CO2 trap CT1 with Vent filter VT1 Kit for 30 / 60 l tank, For LaboStar RO DI, Change interval every 12 months	
LLG06285643	Vent filter for 7 l tank, For LaboStar RO DI, Change interval every 12 months	
LLG06264656	UV lamp 30, 60, 80 liter tank, For Ultra Clear RO EDI, Change interval	1 pz
LLG06285639	Wall bracket ET 60 for 60 l tank, For LaboStar RO DI, Change interval	
LLG06285642	30 l Extension storage tank, For LaboStar RO DI, Change interval	
LLG06285646	Tank 30 l, For All Ultra Clear RO EDI, Change interval	
LLG06285647	Tank 60 l, For All Ultra Clear RO EDI, Change interval	
LLG06285648	Tank 80 l, incl. pump 8 l/min, For Ultra Clear RO EDI 10-30, Change interval	
LLG06285649	Tank 80 l, incl. pump 1.2 l/min, For Ultra Clear RO EDI 10-30, Change interval	
LLG06285651	Tank 80 l, incl. pump 1.2 l/min, For Ultra Clear RO EDI 55, Change interval	
LLG06283462	RO module 55 l/h, For Ultra Clear RO EDI 55 l, Change interval every 2 to 3 years	
LLG06285640	Degassing unit (max. 150 l/h), For LaboStar RO DI, Change interval	
LLG06285650	Tank 80 l, incl. pump 8 l/min, For Ultra Clear RO EDI 55, Change interval	
LLG07981287	Wall bracket ET 30 for 30 l tank, For LaboStar RO DI, Change interval	





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

