

**Cool packs Icecatch****Icecatch****Cool packs Icecatch®**

Icecatch® cool packs for temperature sensitive products used in the pharmaceutical, biotech, medical and food industries. Icecatch® cool packs distinguish themselves through the following specific features:

- Reliable cold storage through the use of immobilised coolant
- For repeat usage due to the use of robust PA/PE film
- Disposal with the domestic trash due to the use of non-poisonous ingredients
- Unobjectionable use in the pharmaceutical and food industries
- Phase transition 0°C

Icecatch®-Gel

- Spreading of heat absorption around the 0°C range through the use of immobilised cooling gel
- Optimal price / heat absorption capacity relationship
- Filled with cooling gel

Icecatch®-Solid

- Maximum dimensional stability during the freezing process through the use of hard foam
- Filled with hard foam

Icecatch®-Solid Insulated 2 - 8 °C

- A perfect cool pack in combination with our System packaging for the +2 °C up to +8 °C transport
- Maximum dimensional stability during the freezing process through the use of hard foam, coolant and integrated isolation layer
- Filled with hard foam and isolation layer



Code	Description	Packaging
LLG06262659	Cool packs Icecatch®, Type Icecatch®-Solid, With hard foam/cool medium, dim. (W x D x H) 195 x 130 x 20 mm, Weight 500 g	27 pz.
LLG06262660	Cool packs Icecatch®, Type Icecatch®-Solid, With hard foam/cool medium, dim. (W x D x H) 105 x 180 x 40 mm, Weight 630 g	12 pz.
LLG09695007	Cool packs Icecatch®, Type Icecatch®-Solid Insulated, With hard foam/cool medium and isolation layer, dim. (W x D x H) 280 x 190 x 30 mm, Weight 1100 g	10 pz
LLG09695000	Cool packs Icecatch®, Type Icecatch®-Gel, With cooling gel, dim. (W x D x H) 90 x 110 x 20 mm, Weight 90 g	120 pz
LLG09695001	Cool packs Icecatch®, Type Icecatch®-Gel, With cooling gel, dim. (W x D x H) 90 x 150 x 20 mm, Weight 170 g	80 pz
LLG09695002	Cool packs Icecatch®, Type Icecatch®-Gel, With cooling gel, dim. (W x D x H) 90 x 180 x 20 mm, Weight 230 g	60 pz
LLG09695004	Cool packs Icecatch®, Type Icecatch®-Gel, With cooling gel, dim. (W x D x H) 140 x 240 x 40 mm, Weight 690 g	24 pz
LLG09695006	Cool packs Icecatch®, Type Icecatch®-Solid, With hard foam/cool medium, dim. (W x D x H) 280 x 190 x 25 mm, Weight 1100 g	10 pz
LLG09695008	Cool packs Icecatch®, Type Icecatch®-Solid Insulated, With hard foam/cool medium and isolation layer, dim. (W x D x H) 195 x 130 x 23 mm, Weight 500 g	24 pz.
LLG09695009	Cool packs Icecatch®, Type Icecatch®-Solid Insulated, With hard foam/cool medium and isolation layer, dim. (W x D x H) 280 x 190 x 45 mm, Weight 1700 g	7 pz.



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

