



Space saving jerrycan Flachmann LaboPlast, PP

Space saving jerrycan Flachmann LaboPlast®, PP

Transparent. With its low width the compact jerrycan Flachman is ideal as a space-saving canister. Different liquids can be stored and decanted in a small space. Flachmann can also be stored lying on its side e.g. in a refrigerator. It is very stable due to ribbed structure and thick sides. The compact jerrycan can also be used as a settling canister. With the optional support stands for 5 l jerrycans, they can be placed safely in the open.

- With or without threaded connector 3/4"
- With tamper-evident cap, PP, blue
- Canister without valve autoclavable
- Graduations: 250ml
- With moulded, indestructible litre graduation

Space-saving Compact set: 1 Flachmann LaboPlast® canister 5 l or 10 l and 1 stopcock Compact (9.140 003)

The Flachmann canister is ideal to store and decant liquids in a small space. The Compact stopcock for canisters and barrels features an impressive compact design, which does not overlap the footprint of the canister and therefore can be stored on its side e.g. within a refrigerator or rack without screwing off the stopcock.



Space-saving StopCock set: 1 Flachmann LaboPlast® canister 5 l or 10 l with screw neck and 1 stopcock (9.139 997)



Code	Description	Packaging
LLG09139996	Space saving jerrycan Flachmann LaboPlast®, PP, l, Support stand for 3 jerrycans (only for 5 l canister), dim. (W x D x H) mm	1 pz
LLG09139999	Space saving jerrycan Flachmann LaboPlast®, PP, l, Support stand for 2 jerrycans (only for 5 l canister), dim. (W x D x H) mm	1 pz
LLG09140001	Space saving jerrycan Flachmann LaboPlast®, PP, 5 l, Without screw neck, dim. (W x D x H) 65 x 335 x 335 mm	1 pz
LLG09140002	Space saving jerrycan Flachmann LaboPlast®, PP, 5 l, With screw neck,, dim. (W x D x H) 65 x 335 x 335 mm	1 pz
LLG09140004	Space saving jerrycan Flachmann LaboPlast®, PP, 5 l, Space-saving Set Compact, dim. (W x D x H) 65 x 335 x 335 mm	1 pz
LLG09140011	Space saving jerrycan Flachmann LaboPlast®, PP, 5 l, Space-saving Set StopCock, dim. (W x D x H) 65 x 335 x 335 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

