

**Tube Rack Rota-Rack, PP****Heathrow Scientific****Tube Rack Rota-Rack®, PP**

Interlocked racks. Rotate each of the four modules independently. Use several sizes of test tubes, culture tubes, centrifuge tubes or microcentrifuge tubes at the same time. Autoclavable. Each module of the Rota-Rack® for large tubes holds 2 x 50ml tubes, 8 x 15ml tubes, 6 x 20mm tubes or 10 x 10/13mm tubes. Each module of the Rota-Rack® for microtubes holds 6 x 15ml tubes, 9 x 1.5/2.0ml tubes, 12 x 0.5ml tubes or 32 x 0.2ml PCR tubes or 4 x 8-tube strips.



Code	Description	Packaging
LLG09193988	Tube Rack Rota-Rack®, PP, Rota-Rack® for large tubes, Width 110 mm, Length 224 mm, Height 110 mm	1 pz.
LLG09193989	Tube Rack Rota-Rack®, PP, Rota-Rack® for micro tubes, Width 100 mm, Length 210 mm, Height 100 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

