



## Industrial jerrycans, HDPE

### Industrial jerrycans, HDPE

Jerrycans with capacity of 2 to 30L.

**Please order closures separately.**



| Code        | Description                                                                                                                              | Packaging |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| LLG09140021 | Industrial jerrycans, HDPE, 2 L, Colour natural, Width 145 mm, Depth 109 mm, Height 196 mm, Weight 110 g, Thread external diam. 45 mm    | 1 pz.     |
| LLG09140030 | Industrial jerrycans, HDPE, 10 L, Colour natural, Width 230 mm, Depth 190 mm, Height 310 mm, Weight 450 g, Thread external diam. 50* mm  | 1 pz.     |
| LLG09140033 | Industrial jerrycans, HDPE, 30 L, Colour blue, Width 320 mm, Depth 284 mm, Height 435 mm, Weight 1250 g, Thread external diam. 60* mm    | 1 pz.     |
| LLG09140038 | Industrial jerrycans, HDPE, 30 L, Colour natural, Width 320 mm, Depth 284 mm, Height 435 mm, Weight 1250 g, Thread external diam. 60* mm | 1 pz.     |
| LLG09140027 | Industrial jerrycans, HDPE, 8 L, Colour natural, Width 230 mm, Depth 190 mm, Height 270 mm, Weight 320 g, Thread external diam. 45 mm    | 1 pz.     |
| LLG09140026 | Industrial jerrycans, HDPE, 6 L, Colour natural, Width 192 mm, Depth 145 mm, Height 294 mm, Weight 265 g, Thread external diam. 45 mm    | 1 pz.     |
| LLG09140022 | Industrial jerrycans, HDPE, 2.5 L, Colour natural, Width 153 mm, Depth 115 mm, Height 208 mm, Weight 125 g, Thread external diam. 45* mm | 1 pz.     |
| LLG09140031 | Industrial jerrycans, HDPE, 25 L, Colour blue, Width 290 mm, Depth 246 mm, Height 450 mm, Weight 1100 g, Thread external diam. 60* mm    | 1 pz.     |
| LLG09140023 | Industrial jerrycans, HDPE, 3 L, Colour natural, Width 153 mm, Depth 115 mm, Height 240 mm, Weight 140 g, Thread external diam. 45* mm   | 1 pz.     |
| LLG09140032 | Industrial jerrycans, HDPE, 20 L, Colour blue, Width 290 mm, Depth 246 mm, Height 385 mm, Weight 950 g, Thread external diam. 60* mm     | 1 pz.     |
| LLG09140037 | Industrial jerrycans, HDPE, 20 L, Colour natural, Width 290 mm, Depth 246 mm, Height 385 mm, Weight 950 g, Thread external diam. 60* mm  | 1 pz.     |
| LLG09140025 | Industrial jerrycans, HDPE, 5 L, Colour natural, Width 190 mm, Depth 145 mm, Height 251 mm, Weight 250 g, Thread external diam. 50* mm   | 1 pz.     |





CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA  
Reagents operates with  
a Certified Quality  
Management System*

