

**Multi-function tumbling shaker PS-M3D****Grant****Multi-function tumbling shaker PS-M3D**

With consistent and reproducible agitation throughout the speed range, multi segment program.

Specifications

Type of movement: orbital (diam. 22mm)

reciprocal (0 to 360°, 30° steps)

vibrating (0 to 5°, 1° steps)

Tumbling angle: 7°

Speed range: 1 to 100rpm

Max. Load: 1.0kg

Timer Type of movement orbital and reciprocal / vibrating : 0 to 250 sec / 0 to 5 sec

Ambient Temperature: +4 to 40°C

Dimensions Plattform (W x D): 215 x 215 mm

Dimensions (W x D x H): 235 x 235 x 140 mm

Weight: 2kg

Power supply: 12V



Code	Description	Packaging
LLG09721017	Multi-function rotator 3D mixer incl.EU-/UK-plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

