

**Vortex mixer PV-1****Grant****Vortex mixer PV-1**

With low profile and small footprint for gentle mixing to vigorous re-suspension of cell or chemical pellets.

- variable speed: 750rpm to 3000rpm
- continuous or touch operation
- for tubes up to 28.5 mm (50 ml) diameter

Specifications

Type of movement: orbital

Orbit: 4mm

Speed range: 750 to 3000rpm

Dimensions (WxDxH): 90 x 150 x 80mm

Temperature range: 4 to 40°C



Code	Description	Packaging
LLG09721000	Mini-Vortex mixer PV-1 for tubes 1,5 ml - 50 ml incl. EU-/UK-plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

