

**Centrifuge Multi-spin PCV-6000****Grant****Centrifuge Multi-spin PCV-6000**

Combined centrifuge/vortex mixer for microtubes and 0.2ml strips. For independent centrifuge and vortex mixing operations, or a programmable sequence of spin and mix phases of multiple tubes.

- Vortex mixing intensity: three stage adjustment (soft, medium and hard)
- Interchangeable rotors for 2.0ml, 1.5ml, 0.5ml, 0.2ml microtube combinations and 0.2ml strips
- Timed vortex mixing: 1sec to 20sec
- Supplied with 2 rotors: 12 x 1.5/2.0ml and 12 x 0.5ml plus 12 x 0.2ml

Specifications

Max. RCF(bottom of tubes) 2350xg

Max. speed 6000rpm

Max. capacity: 12 x 1.5/2.0ml or 12 x 0.5ml plus 12 x 0.2ml or 2 x 0.2ml strips of 8

Spin Timer 1sec to 30min

Power consumption 24W

Dimensions (W x D x H): 190 x 235 x 125mm

Weight: 2.5kg

Power supply: 100 - 240V, 50/60Hz

Input voltage 24V



Code	Description	Packaging
LLG09721004	Centrifuge Multi-spin PCV-6000, Type PCV-6000, Centrifuge Multi-spin	1 pz.
LLG09721008	Accessories for Centrifuges PCV-2400 and PCV-6000, Type PSR-16, Rotor for 2 x 8-well 0.2 ml microtube strips	1 pz.
LLG09721007	Accessories for Centrifuges PCV-2400 and PCV-6000, Type PR2-05-02, Rotor for 6 x 1.5/2.0 ml + 6 x 0.2 ml tubes	1 pz.
LLG09721006	Accessories for Centrifuges PCV-2400 and PCV-6000, Type PR2-05, Rotor for 8 x 1.5/2.0 ml + 8 x 0.5 ml tubes	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

