



Cryogenic boxes, 81 well, PP, autoclavable

Cryogenic boxes, 81 well, PP, autoclavable

The boxes are made of polypropylene and can be used in temperatures as low as -90 °C. The boxes are autoclavable at +121 °C for 20 minutes. The boxes have a fixed grid divider with 9x9 cells and are suited for 1.0 ml and 2 ml cryogenic vials. Cellcodes are printed on the lid for easy identification.



Code	Description	Packaging
LLG09698708	Cryogenic boxes, 81 well, PP, autoclavable, Type Blue, dim. (W x D x H) 130 x 130 x 50 mm	1 pz.
LLG09698709	Cryogenic boxes, 81 well, PP, autoclavable, Type Green, dim. (W x D x H) 130 x 130 x 50 mm	1 pz.
LLG09698710	Cryogenic boxes, 81 well, PP, autoclavable, Type Natural, dim. (W x D x H) 130 x 130 x 50 mm	1 pz.
LLG09698711	Cryogenic boxes, 81 well, PP, autoclavable, Type Red, dim. (W x D x H) 130 x 130 x 50 mm	1 pz.
LLG09698712	Cryogenic boxes, 81 well, PP, autoclavable, Type Yellow, dim. (W x D x H) 130 x 130 x 50 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

