

**Soxhlet extraction stillheads, DURAN tubing****Lenz****Soxhlet extraction stillheads, DURAN® tubing**

Soxhlet extraction stillheads, made of DURAN® tubing. Extractor heads acc. to DIN 12602 and in addition to DIN standard (*).



Code	Description	Packaging
LLG09043018	Soxhlet extraction stillheads, DURAN® tubing, For extractor 200* ml, Cone 29/32 NS, Socket 45/40 NS	1 pz.
LLG09043015	Soxhlet extraction stillheads, DURAN® tubing, For extractor 250 ml, Cone 29/32 NS, Socket 45/40 NS	1 pz.
LLG09043012	Soxhlet extraction stillheads, DURAN® tubing, For extractor 70* ml, Cone 29/32 NS, Socket 34/35 NS	1 pz.
LLG09043011	Soxhlet extraction stillheads, DURAN® tubing, For extractor 30 ml, Cone 29/32 NS, Socket 29/32 NS	1 pz.
LLG09043013	Soxhlet extraction stillheads, DURAN® tubing, For extractor 100 ml, Cone 29/32 NS, Socket 45/40 NS	1 pz.
LLG09043014	Soxhlet extraction stillheads, DURAN® tubing, For extractor 150* ml, Cone 29/32 NS, Socket 45/40 NS	1 pz.
LLG09043016	Soxhlet extraction stillheads, DURAN® tubing, For extractor 500 ml, Cone 29/32 NS, Socket 60/46 NS	1 pz.
LLG09043017	Soxhlet extraction stillheads, DURAN® tubing, For extractor 1000 ml, Cone 29/32 NS, Socket 71/51 NS	1 pz.
LLG09043019	Soxhlet extraction stillheads, DURAN® tubing, For extractor 300* ml, Cone 29/32 NS, Socket 60/46 NS	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

