

**Cryogenic boxes, PP, 133 x 133****Ratiolab****Cryogenic boxes, PP, 133 x 133**

For cryogenic storage of cryotubes.

- 9 x 9 grid
- Grid numerically coded
- Removable grid inserts
- Hinge with safe snap-on lid
- Stackable
- Drain bores at the bottom for dew liquid
- Temperature resistant from -90 to 121°C
- Autoclavable



| Code        | Description                                                                          | Packaging |
|-------------|--------------------------------------------------------------------------------------|-----------|
| LLG09405900 | Cryogenic boxes, PP, 133 x 133, For tubes 1.2 ... 2 ml, Colour Black, Height 52 mm   | 1 pz.     |
| LLG09405901 | Cryogenic boxes, PP, 133 x 133, For tubes 3 ... 4 ml, Colour Black, Height 75 mm     | 1 pz.     |
| LLG09405898 | Cryogenic boxes, PP, 133 x 133, For tubes 1.2 ... 2 ml, Colour Blue, Height 52 mm    | 1 pz.     |
| LLG09405899 | Cryogenic boxes, PP, 133 x 133, For tubes 3 ... 4 ml, Colour Blue, Height 75 mm      | 1 pz.     |
| LLG09405896 | Cryogenic boxes, PP, 133 x 133, For tubes 1.2 ... 2 ml, Colour Green, Height 52 mm   | 1 pz.     |
| LLG09405897 | Cryogenic boxes, PP, 133 x 133, For tubes 3 ... 4 ml, Colour Green, Height 75 mm     | 1 pz.     |
| LLG09405890 | Cryogenic boxes, PP, 133 x 133, For tubes 1.2 ... 2 ml, Colour Natural, Height 52 mm | 1 pz.     |
| LLG09405891 | Cryogenic boxes, PP, 133 x 133, For tubes 3 ... 4 ml, Colour Natural, Height 75 mm   | 1 pz.     |
| LLG09405894 | Cryogenic boxes, PP, 133 x 133, For tubes 1.2 ... 2 ml, Colour Red, Height 52 mm     | 1 pz.     |
| LLG09405895 | Cryogenic boxes, PP, 133 x 133, For tubes 3 ... 4 ml, Colour Red, Height 75 mm       | 1 pz.     |
| LLG09405892 | Cryogenic boxes, PP, 133 x 133, For tubes 1.2 ... 2 ml, Colour Yellow, Height 52 mm  | 1 pz.     |
| LLG09405893 | Cryogenic boxes, PP, 133 x 133, For tubes 3 ... 4 ml, Colour Yellow, Height 75 mm    | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

