

**Flexible Distributors, PP****BOLA****Flexible Distributors, PP**

For connecting elastic tubing by simply slipping it onto hose connectors. Suitable for GL45 threaded bottles. With screw cap made of glass-fibre reinforced PP, ideal for transferring liquids or gases. The angled tubing connectors can be attached to a liquid source, a pressure source or a vacuum system. Avoid crimping of the tubing. The distributor can be loosened by simply turning the screw cap. The complete distribution system can then be fixed onto another bottle.



Code	Description	Packaging
LLG07634892	Flexible Distributors, PP, For tubing bore 6 - 9 mm, number of nozzles 2	1 pz.
LLG09209858	Flexible Distributors, PP, For tubing bore 6 - 9 mm, number of nozzles 3	1 pz.
LLG09209859	Flexible Distributors, PP, For tubing bore 7 - 11 mm, number of nozzles 2	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

