



Disposable samplers for solids, HDPE

Bürkle

Disposable samplers for solids, HDPE

In numerous situations, rapid and contamination-free samples have to be taken. DispoSamplers are perfect for this application. They are cleanroom manufactured, individually packed and always ready to use. Avoid expensive decontamination or cleaning of soiled samplers. Gamma sterilised versions are also available.

MultiDispo: For free-flowing powders and small-grain granulates. Easy emptying through opening in the handle.

ZoneDispo: For free-flowing powders and small-grain granulates. Collected sample can be transported in closed sampler.

TargetDispo: For taking target-point samples from powders and granulates.

PowderDispo: For taking cross-section samples of powdery or slightly cohesive substances.

DispoPicker: For taking cross-section samples from powders and small-grain granulates. With closed tip.

DispoLance: For direct removal and filling of free-flowing powders and small-grain granulates.

MicroDispo: For free-flowing powders and small-grain bulk materials. For small sample quantities (10 ml).



Code	Description	Packaging
LLG09303728	Disposable samplers for solids, HDPE, Type MicroDispo, non-sterile, 10 ml, Length 500 mm	20 pz.
LLG09303729	Disposable samplers for solids, HDPE, Type MicroDispo, sterile, 10 ml, Length 500 mm	20 pz
LLG09303730	Disposable samplers for solids, HDPE, Type MicroDispo, non-sterile, 10 ml, Length 1000 mm	20 pz.
LLG09303731	Disposable samplers for solids, HDPE, Type MicroDispo, sterile, 10 ml, Length 1000 mm	20 pz
LLG06259428	Disposable samplers for solids, HDPE, Type PowderDispo, sterile, 100 ml, Length 500 mm	20 pz
LLG06269015	Disposable samplers for solids, HDPE, Type ZoneDispo, sterile, 100 ml, Length 500 mm	20 pz
LLG06282439	Disposable samplers for solids, HDPE, Type MultiDispo, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG06282440	Disposable samplers for solids, HDPE, Type MultiDispo, sterile, 100 ml, Length 500 mm	20 pz
LLG06282444	Disposable samplers for solids, HDPE, Type TargetDispo, sterile, 100 ml, Length 500 mm	20 pz
LLG06282456	Disposable samplers for solids, HDPE, Type DispoPicker, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG06282457	Disposable samplers for solids, HDPE, Type DispoPicker, sterile, 100 ml, Length 500 mm	20 pz
LLG07655427	Disposable samplers for solids, HDPE, Type ZoneDispo, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG09303720	Disposable liquid samplers, HDPE, Type LiquiDispo, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG09303732	Disposable liquid samplers, HDPE, Type ViscoDispo, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG09303736	Disposable samplers for solids, HDPE, Type PowderDispo, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG06282438	Disposable samplers for solids, HDPE, Type ZoneDispo, sterile, 200 ml, Length 1000 mm	20 pz
LLG06282441	Disposable samplers for solids, HDPE, Type MultiDispo, non-sterile, 200 ml, Length 1000 mm	20 pz.
LLG06282442	Disposable samplers for solids, HDPE, Type MultiDispo, sterile, 200 ml, Length 1000 mm	20 pz
LLG06282446	Disposable samplers for solids, HDPE, Type TargetDispo, sterile, 200 ml, Length 1000 mm	20 pz
LLG06282449	Disposable liquid samplers, HDPE, Type DispoTube, non-sterile, 200 ml, Length 1000 mm	20 pz



Code	Description	Packaging
LLG09303722	Disposable liquid samplers, HDPE, Type LiquiDispo, non-sterile, 200 ml, Length 1000 mm	20 pz.
LLG09303734	Disposable liquid samplers, HDPE, Type ViscoDispo, non-sterile, 200 ml, Length 1000 mm	20 pz.
LLG06282452	Disposable samplers for solids, HDPE, Type PowderDispo, non-sterile, 60 ml, Length 300 mm	20 pz.
LLG06282453	Disposable samplers for solids, HDPE, Type PowderDispo, sterile, 60 ml, Length 300 mm	20 pz
LLG06282454	Disposable samplers for solids, HDPE, Type DispoPicker, non-sterile, 60 ml, Length 300 mm	20 pz.
LLG06282455	Disposable samplers for solids, HDPE, Type DispoPicker, sterile, 60 ml, Length 300 mm	20 pz
LLG06282458	Disposable samplers for solids, HDPE, Type DispoLance, non-sterile, ml, Length 300 mm	20 pz.
LLG06282459	Disposable samplers for solids, HDPE, Type DispoLance, sterile, ml, Length 300 mm	20 pz
LLG06282460	Disposable samplers for solids, HDPE, Type DispoLance, non-sterile, ml, Length 500 mm	20 pz.
LLG06282461	Disposable samplers for solids, HDPE, Type DispoLance, sterile, ml, Length 500 mm	20 pz
LLG06282443	Disposable samplers for solids, HDPE, Type TargetDispo, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG06282445	Disposable samplers for solids, HDPE, Type TargetDispo, non-sterile, 200 ml, Length 1000 mm	20 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

