

**Magnetic stirrer with heating LLG-uniSTIRRER 3 pro, complete set****LLG Labware****Magnetic stirrer with heating LLG-uniSTIRRER 3 pro, complete set**

Universal magnetic stirrer with outstanding price-performance ratio and new, improved design.

- High level of safety thanks to **"Hot" display** with continuing stirring function until the plate temperature has dropped below 50 °C, even when the device is switched off.
- Integrated control function for precise control of the media temperature by a connected PT1000 temperature probe.
- Stainless steel plate with ceramic coating
- Digital speed and temperature control
- LED display
- Fixed safety circuit of the hotplate temperature at 320 °C

**Scope of supply:** Magnetic stirrer, PT1000 temperature probe, rod, holder

**Specifications :**

- Stirring quantity: 3 l
- Speed range : 200...1500 rpm
- Display accuracy :  $\pm 10$  rpm
- Magnetic bar : max. 50 mm
- Display: LED
- Heating plate :  $\varnothing$  135 mm
- Permissible ambient temp. : 5 ... 40 °C
- Temperature range : RT ...+280 °C
- Setting accuracy :  $\pm 1$  °C
- Control accuracy :  $\pm 2$  °C
- Heat output : 500 W
- Safety temperature : 320 °C
- Dimensions (W x D x H): 155x215x100mm
- Weight: 1,5 kg
- Power supply : 200...240V, 50/60Hz
- IP code: IP 21
- warranty 3 years



Code	Description	Packaging
LLG06263400	LLG-uniSTIRRER 3 pro heating stirrer + PT1000 / EU plug	1 pz
LLG06263401	LLG-uniSTIRRER 3 pro heating stirrer + PT1000 / UK plug	



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

