



Grinding bowls and grinding balls for mini mill PULVERISETTE 23

Fritsch

Grinding bowls and grinding balls for mini mill PULVERISETTE 23

Grinding bowl|5ml|10ml|15ml

Useful capacity (sample volume)|0.1 - 1ml|0.2 - 1ml|0.5 - 5ml

Ball diameter|||

5mm|20|30|60

10mm|1|3|5

15mm|-|1|2



Code	Description	Packaging
LLG09739164	Grinding bowls and grinding balls for mini mill PULVERISETTE 23, Type Grinding ball, 15 mm, hardened, stainless steel	
LLG09739349	Grinding bowls and grinding balls for mini mill PULVERISETTE 23, Type Grinding bowl, hardened, stainless steel, 10 ml	1 pz.
LLG09739348	Grinding bowls and grinding balls for mini mill PULVERISETTE 23, Type Grinding bowl, hardened, stainless steel, 15ml	1 pz.
LLG09739343	Grinding bowls and grinding balls for mini mill PULVERISETTE 23, Type Grinding bowl, PTFE, 5ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

