



Accessories for Ultrasonic Processor UP50H/UP100H


Accessories for Ultrasonic Processor UP50H/UP100H

MS0.5 - MS10: Sonotrodes made of titanium, approx. length 80 mm (MS10 only for UP100H).

SB1 + SB2: Sound protection box with height adjustable desk top and support rod dia. 16 mm, for UP200S or UP400S (SB2-16 only for UP50H or UP100H).

ST 1: Stand, pole diam. 16 mm made of stainless steel, base L x W x H: 300 x 150 x 600 mm.

T1: Timer from 00:00 to 99:59 (min:sec).



Code	Description	Packaging
LLG09106602	Accessories for Ultrasonic Processor UP50H/UP100H, Type MS1, Tip Ø 1 mm, For volume 0,1 - 5,0 ml	1 pz.
LLG09106604	Accessories for Ultrasonic Processor UP50H/UP100H, Type MS7, Tip Ø 7 mm, For volume 10 - 250 ml	1 pz.
LLG04005990	Accessories for Ultrasonic Processor UP50H/UP100H, Type MS10, Tip Ø 10 mm, For volume 20 - 500 ml	1 pz.
LLG04005989	Accessories for Ultrasonic Processor UP50H/UP100H, Type MS3, Tip Ø 3 mm, For volume 5 - 100 ml	1 pz.
LLG06226573	Accessories for Ultrasonic Processor UP50H/UP100H, Type ST 1, Tip -, For volume - ml	1 pz.
LLG06235003	Accessories for Ultrasonic Processor UP50H/UP100H, Type T1, Tip -, For volume - ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

