



Sonotrode for compact ultrasonic processor UP200S

Hielscher Ultrasonics

Sonotrode for compact ultrasonic processor UP200S



Code	Description	Packaging
LLG09106607	Sonotrode for compact ultrasonic processor UP200S, Type S1, Tip dia. 1 mm, For volume 0.1 ... 5 ml ml	1 pz.
LLG09106609	Sonotrode for compact ultrasonic processor UP200S, Type S7, Tip dia. 7 mm, For volume 20 ... 500 ml ml	1 pz.
LLG09106610	Sonotrode for compact ultrasonic processor UP200S, Type S14, Tip dia. 14 mm, For volume 50 ... 1000 ml ml	1 pz.
LLG09106608	Sonotrode for compact ultrasonic processor UP200S, Type S3, Tip dia. 3 mm, For volume 5 ... 200 ml ml	1 pz.
LLG09106611	Sonotrode for compact ultrasonic processor UP200S, Type S14D*, Tip, For volume 5 ... 25 l/h ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

