

**Dispersing heads for Ultra Turrax T18 digital****IKA****Dispersing heads for Ultra Turrax® T18 digital**

Plastic dispersing heads: Autoclavable

Stainless steel dispersing heads: Ultimate fineness emulsions 1 to 10µm and suspensions 10 to 50µm, sterilisable by all methods.



| Code        | Description                                                                                                                            | Packaging |
|-------------|----------------------------------------------------------------------------------------------------------------------------------------|-----------|
| LLG09570842 | Dispersing heads for Ultra Turrax® T18 digital, Type S 18 N - 19 G,<br>Working volume 10 to 1500 ml, Diam. rotor / stator 12.7 / 19 mm | 1 pz.     |
| LLG09570841 | Dispersing heads for Ultra Turrax® T18 digital, Type S 18 N - 10 G,<br>Working volume 1 to 100 ml, Diam. rotor / stator 7.5 / 10 mm    | 1 pz.     |
| LLG09570857 | Dispersing heads for Ultra Turrax® T18 digital, Type S 18 D-14 G-KS,<br>Working volume 10 to 500 ml, Diam. rotor / stator 9.5 / 14 mm  | 10 pz.    |
| LLG09570855 | Dispersing heads for Ultra Turrax® T18 digital, Type S 18 D-10 G-KS,<br>Working volume 10 to 100 ml, Diam. rotor / stator 6.75 / 10 mm | 10 pz.    |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

