

**Dispersing heads for Ultra Turrax T18 digital****IKA****Dispersing heads for Ultra Turrax® T18 digital**

Plastic dispersing heads: Autoclavable

Stainless steel dispersing heads: Ultimate fineness emulsions 1 to 10µm and suspensions 10 to 50µm, sterilisable by all methods.



Code	Description	Packaging
LLG09570842	Dispersing heads for Ultra Turrax® T18 digital, Type S 18 N - 19 G, Working volume 10 to 1500 ml, Diam. rotor / stator 12.7 / 19 mm	1 pz.
LLG09570841	Dispersing heads for Ultra Turrax® T18 digital, Type S 18 N - 10 G, Working volume 1 to 100 ml, Diam. rotor / stator 7.5 / 10 mm	1 pz.
LLG09570857	Dispersing heads for Ultra Turrax® T18 digital, Type S 18 D-14 G-KS, Working volume 10 to 500 ml, Diam. rotor / stator 9.5 / 14 mm	10 pz.
LLG09570855	Dispersing heads for Ultra Turrax® T18 digital, Type S 18 D-10 G-KS, Working volume 10 to 100 ml, Diam. rotor / stator 6.75 / 10 mm	10 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

