



Accessories for Grinding mill, A11 basic

IKA

Accessories for Grinding mill, A11 basic

A 11.1 Spare beater: For pulverising substances with a Mohs hardness up to 6.

A 11.2 Cutter: For pulverising soft, fibrous grinding materials.

A 11.3 Beater, hard metal: For pulverising substances with a Mohs hardness up to 9, chrome-carbide coated

A 11.4 Grinding container, 250ml: Made of polycarbonate with stainless steel inlet. Not suitable for cooling with N₂, only applicable with double beater A 11.6.

A 11.5 Spare grinding container, 80ml: Made of Tefcel (PTFE, glass fiber-reinforced) with stainless steel inlet. Excellent resistance to chemicals and low temperatures (- 200 °C).

A 11.6 Dubbel beater: For use up to Mohs hardness 3. Only applicable with grinding chamber A 11.4.

A 11.7 Funnel: Prevents splashing by pouring in liquid nitrogen in the grinding chamber A 11.5.



Code	Description	Packaging
LLG09737807	Accessories for Grinding mill, A11 basic, Type A 11.3, Beater	1 pz.
LLG09737817	Accessories for Grinding mill, A11 basic, Type A 11.6, Double beater	1 pz.
LLG09737806	Accessories for Grinding mill, A11 basic, Type A 11.2, Cutter	1 pz.
LLG09737818	Accessories for Grinding mill, A11 basic, Type A 11.7, Funnel	1 pz.
LLG09737805	Accessories for Grinding mill, A11 basic, Type A 11.4, Grinding container, 250ml	1 pz.
LLG09737815	Accessories for Grinding mill, A11 basic, Type A 11.1, Spare beater	1 pz.
LLG09737816	Accessories for Grinding mill, A11 basic, Type A 11.5, Spare grinding container, 80 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

