

**Erlenmeyer flasks, DURAN, with screw neck****DWK Life Sciences****Erlenmeyer flasks, DURAN®, with screw neck**

With DIN screw neck, PBT-screwthread cap (autoclavable up to 180 °C), and approximate graduations. The flask is suitable for storage and cultivation. Glass Type I / neutral glass as per USP, EP and JP. Autoclavable; undo cap before autoclaving. With retrace code (Batch Identification). Certificate via internet.



Code	Description	Packaging
LLG09141224	Erlenmeyer flasks, DURAN®, with screw neck, 100 ml, Bottom diam. 64 mm, Height 109 mm, Neck thread 25 GL	1 pz
LLG09141244	Erlenmeyer flasks, DURAN®, with screw neck, 500 ml, Bottom diam. 105 mm, Height 180 mm, Neck thread 32 GL	1 pz
LLG09141254	Erlenmeyer flasks, DURAN®, with screw neck, 1000 ml, Bottom diam. 131 mm, Height 225 mm, Neck thread 32 GL	1 pz
LLG09141236	Erlenmeyer flasks, DURAN®, with screw neck, 250 ml, Bottom diam. 85 mm, Height 149 mm, Neck thread 32 GL	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

