

**Gel electrophoresis tank OmniPage Mini****Cleaver Scientific****Gel electrophoresis tank OmniPage Mini**

Vertical chamber especially for protein analysis in up to 4 gels. The special manufacturing process results in a device with unique quality and durability at a very reasonable price. Easy handling and fast cooling. The dimensions are compatible with all 8 x 10 and 10 x 10 cm ready-to-use gels. High sample throughput by up to 4 gels with 20 samples each. Low buffer volume from 250 to 1200 ml.

Supplied with:

OmniPage Mini: chamber with platinum electrodes, 2 sets of glass plates with 1 mm thick bonded Spacers, electrodes, 2 x 12 sample combs.

OmniPage Mini SYS: OmniPage Mini + caster.

OmniPage Mini CBS: OmniPage Mini + caster and blotting insert.

Technical specifications

Gel dimensions (W x L): 7.5 x 8 cm

Unit dimensions (W x D x H): 19 x 13 x 15 cm

Max. sample capacity: 80 samples, 20 samples per gel

Buffer volume: min. 250 ml; max 1200 ml



Code	Description	Packaging
LLG06232290	Accessories for Electrophoresis Tank OmniPage Mini, Comb 10 sample, Thickness 1.00 mm, Sample volume 40 μ l	1 pz.
LLG06232176	Accessories for Electrophoresis Tank OmniPage Mini, Comb 12 sample, Thickness 1.00 mm, Sample volume 35 μ l	1 pz.
LLG09584732	Accessories for Electrophoresis Tank OmniPage Mini, Comb 12 sample, Thickness 0.75 mm, Sample volume 25 μ l	1 pz.
LLG09584734	Accessories for Electrophoresis Tank OmniPage Mini, Comb 12 sample, Thickness 1.50 mm, Sample volume 50 μ l	1 pz.
LLG06232291	Accessories for Electrophoresis Tank OmniPage Mini, Comb 20 sample, Thickness 1.00 mm, Sample volume 20 μ l	1 pz.
LLG06232797	Accessories for Electrophoresis Tank OmniPage Mini, Comb 8 sample MC, Thickness 1.00 mm, Sample volume 60 μ l	1 pz.
LLG09584652	Gel electrophoresis tank OmniPage Mini, OmniPage Mini	1 pz.
LLG06229038	Gel electrophoresis tank OmniPage Mini, OmniPage Mini CBS	1 pz.
LLG09584654	Gel electrophoresis tank OmniPage Mini, OmniPage Mini SYS	1 pz.
LLG09584730	Accessories for Electrophoresis Tank OmniPage Mini, Plain glass plates with bonded spacers, Thickness 0.75 mm, Sample volume μ l	1 pz.
LLG09584731	Accessories for Electrophoresis Tank OmniPage Mini, Plain glass plates with bonded spacers, Thickness 1.50 mm, Sample volume μ l	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
Reagents operates with
a Certified Quality
Management System

