

**Gas washing bottles acc. to Drechsel, complete, DURAN tubing****Gas washing bottles acc. to Drechsel, complete, DURAN® tubing**

Bottle, DURAN®. Head made from DURAN® tubing, with screw connections. Head available with or without porosity 1 filter disc, insert height adjustable, capacity: 500ml, central screw cap: GL 45, lateral screw caps: GL 14.



Code	Description	Packaging
LLG06073489	Gas Washbottle without filter disc, 500 ml	1 pz.
LLG06206086	Gas Washbottle Head with filter disc, Porosity 1	
LLG07635908	Gas Washbottle with filter disc, 500 ml, Porosity 1	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

