

**Nitrogen round-bottom flasks (Schlenk-flasks), DURAN****Lenz****Nitrogen round-bottom flasks (Schlenk-flasks), DURAN®**

With an angled side arm with 2.5 mm glass stopcock. Material: DURAN®



Code	Description	Packaging
LLG06073481	Nitrogen round-bottom flasks (Schlenk-flasks), DURAN®, 250 ml, Socket 29/32 NS, Bore 2.5 mm	1 pz.
LLG06206524	Nitrogen round-bottom flasks (Schlenk-flasks), DURAN®, 250 ml, Socket 14/23 NS, Bore 2.5 mm	1 pz.
LLG06307327	Nitrogen round-bottom flasks (Schlenk-flasks), DURAN®, 100 ml, Socket 14/23 NS, Bore 2.5 mm	1 pz.
LLG06307820	Nitrogen round-bottom flasks (Schlenk-flasks), DURAN®, 50 ml, Socket 14/23 NS, Bore 2.5 mm	1 pz.
LLG06307821	Nitrogen round-bottom flasks (Schlenk-flasks), DURAN®, 100 ml, Socket 29/32 NS, Bore 2.5 mm	1 pz.
LLG06307822	Nitrogen round-bottom flasks (Schlenk-flasks), DURAN®, 500 ml, Socket 29/32 NS, Bore 2.5 mm	1 pz.
LLG06800337	Nitrogen round-bottom flasks (Schlenk-flasks), DURAN®, 25 ml, Socket 14/23 NS, Bore 2.5 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

