

**Extractor heads for specific light solvents, DURAN tubing****Lenz****Extractor heads for specific light solvents, DURAN® tubing**

Extractor heads for the extraction of liquids with specific light solvents, with insert filter (porosity: 0), made of DURAN® tubing.

Distribution tubes with frit (porosity: 0), for specific light solvents, made of DURAN® tubing.



Code	Description	Packaging
LLG06258695	Extractor heads for specific light solvents, DURAN® tubing, Type Extractor head, Extractor 250 ml, Condenser 45/40 NS	1 pz.
LLG06801962	Extractor heads for specific light solvents, DURAN® tubing, Type Extractor head, Extractor 100 ml, Condenser 45/40 NS	1 pz.
LLG06902445	Extractor heads for specific light solvents, DURAN® tubing, Type Extractor head, Extractor 500 ml, Condenser 60/46 NS	1 pz.
LLG06258692	Extractor heads for specific light solvents, DURAN® tubing, Type Distribution tube, Extractor 100 ml, Condenser - NS	1 pz.
LLG06258693	Extractor heads for specific light solvents, DURAN® tubing, Type Distribution tube, Extractor 250 ml, Condenser - NS	1 pz.
LLG06902444	Extractor heads for specific light solvents, DURAN® tubing, Type Distribution tube, Extractor 500 ml, Condenser - NS	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

