



## Condensers acc. to Davies, DURAN tubing

### Condensers acc. to Davies, DURAN® tubing

Made of DURAN® tubing. For condensing vapours, with NS socket and cone, with double sized cooling area, water connection: connection with glass olives or GL 14 threads with two threaded oplastic olives.



Code	Description	Packaging
LLG06228159	Condensers acc. to Davies, DURAN® tubing, Effective length 250 mm, Cone 24/29 NS, Socket 24/29 NS, Connection nozzle, Cooling surface 460 cm <sup>2</sup>	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

