



Condensers acc. to Davies, DURAN tubing

Condensers acc. to Davies, DURAN® tubing

Made of DURAN® tubing. For condensing vapours, with NS socket and cone, with double sized cooling area, water connection: connection with glass olives or GL 14 threads with two threaded plastic olives.



Code	Description	Packaging
LLG06228159	Condensers acc. to Davies, DURAN® tubing, Effective length 250 mm, Cone 24/29 NS, Socket 24/29 NS, Connection nozzle, Cooling surface 460 cm ²	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

