



Silicone rubber stoppers

Silicone rubber stoppers

Transparent. To DIN 12871. Solid. Operating temperature up to 220 °C.
Supplied singly.



Code	Description	Packaging
LLG09230503	Silicone rubber stoppers, Bottom diam. 3.5 mm, Top diam. 6.5 mm, Height 15 mm	1 pz.
LLG09230505	Silicone rubber stoppers, Bottom diam. 5 mm, Top diam. 9 mm, Height 20 mm	1 pz.
LLG09230508	Silicone rubber stoppers, Bottom diam. 8 mm, Top diam. 12 mm, Height 20 mm	1 pz.
LLG09230510	Silicone rubber stoppers, Bottom diam. 10.5 mm, Top diam. 14.5 mm, Height 20 mm	1 pz.
LLG09230512	Silicone rubber stoppers, Bottom diam. 12.5 mm, Top diam. 16.5 mm, Height 20 mm	1 pz.
LLG09230514	Silicone rubber stoppers, Bottom diam. 14 mm, Top diam. 18 mm, Height 20 mm	1 pz.
LLG09230517	Silicone rubber stoppers, Bottom diam. 17 mm, Top diam. 22 mm, Height 25 mm	1 pz.
LLG09230518	Silicone rubber stoppers, Bottom diam. 18 mm, Top diam. 24 mm, Height 30 mm	1 pz.
LLG09230521	Silicone rubber stoppers, Bottom diam. 21 mm, Top diam. 27 mm, Height 30 mm	1 pz.
LLG09230523	Silicone rubber stoppers, Bottom diam. 23 mm, Top diam. 29 mm, Height 30 mm	1 pz.
LLG09230526	Silicone rubber stoppers, Bottom diam. 26 mm, Top diam. 32 mm, Height 30 mm	1 pz.
LLG09230531	Silicone rubber stoppers, Bottom diam. 31 mm, Top diam. 38 mm, Height 35 mm	1 pz.
LLG09230536	Silicone rubber stoppers, Bottom diam. 36 mm, Top diam. 44 mm, Height 40 mm	1 pz.
LLG09230541	Silicone rubber stoppers, Bottom diam. 41 mm, Top diam. 49 mm, Height 40 mm	1 pz.
LLG09230529	Silicone rubber stoppers, Bottom diam. 29 mm, Top diam. 35 mm, Height 30 mm	1 pz.
LLG09230547	Silicone rubber stoppers, Bottom diam. 47 mm, Top diam. 55 mm, Height 40 mm	1 pz.
LLG09230556	Silicone rubber stoppers, Bottom diam. 56 mm, Top diam. 65 mm, Height 45 mm	1 pz.
LLG09230560	Silicone rubber stoppers, Bottom diam. 60 mm, Top diam. 70 mm, Height 50 mm	1 pz.
LLG09230564	Silicone rubber stoppers, Bottom diam. 64.5 mm, Top diam. 75.5 mm, Height 55 mm	1 pz.
LLG09230571	Silicone rubber stoppers, Bottom diam. 71 mm, Top diam. 83 mm, Height 60 mm	1 pz.

Code	Description	Packaging
LLG09230550	Silicone rubber stoppers, Bottom diam. 50.5 mm, Top diam. 59.5 mm, Height 45 mm	1 pz.
LLG09230579	Silicone rubber stoppers, Bottom diam. 79 mm, Top diam. 92 mm, Height 65 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

