

**Mens laboratory coats Type 81996, 100% cotton****Mens laboratory coats Type 81996, 100% cotton**

Reversible collar, concealed buttons, 1 outside breast pocket left, 1 inside breast pocket right, 2 side pockets with diagonal access with light grey piping, 2 short side slits, straight, untapered back, straight sleeve seam, Öko-Tex Standard 100

Material: 100% cotton

Colour: White

Weight: approx. 240g/m²

Length: approx. 96cm



Code	Description	Packaging
LLG09414371	Mens laboratory coats Type 81996, 100% cotton, Size 44/46	1 pz.
LLG09414373	Mens laboratory coats Type 81996, 100% cotton, Size 48/50	1 pz.
LLG09414375	Mens laboratory coats Type 81996, 100% cotton, Size 52/54	1 pz.
LLG09414377	Mens laboratory coats Type 81996, 100% cotton, Size 56/58	1 pz.
LLG09414379	Mens laboratory coats Type 81996, 100% cotton, Size 60/62	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

