

**Reaction flasks, flat flange, DURAN****DWK Life Sciences****Reaction flasks, flat flange, DURAN®**

Flanged neck with groove for vacuum use. 500ml and 1000ml vessels are completely cylindrical. Capacities indicated are brim-full volumes. Glass Type I / neutral glass as per USP, EP and JP. Autoclavable.



Code	Description	Packaging
LLG09142024	Reaction flasks, flat flange, DURAN®, 100 ml, 195 ml, Flange diam. 100 mm, Body diam. 70 mm, Height 85 mm, Pressure 2.5* bar, DN 60	1 pz
LLG09142036	Reaction flasks, flat flange, DURAN®, 250 ml, 315 ml, Flange diam. 100 mm, Body diam. 70 mm, Height 125 mm, Pressure 2.5* bar, DN 60	1 pz
LLG09142044	Reaction flasks, flat flange, DURAN®, 500 ml, 740 ml, Flange diam. 138 mm, Body diam. 106 mm, Height 120 mm, Pressure 1.5* bar, DN 100	1 pz
LLG09142054	Reaction flasks, flat flange, DURAN®, 1000 ml, 1395 ml, Flange diam. 138 mm, Body diam. 106 mm, Height 205 mm, Pressure 1.5* bar, DN 100	1 pz
LLG09142071	Reaction flasks, flat flange, DURAN®, 4000 ml, 5765 ml, Flange diam. 184 mm, Body diam. 200 mm, Height 290 mm, Pressure 1.0* bar, DN 150	1 pz
LLG09142063	Reaction flasks, flat flange, DURAN®, 2000 ml, 2620 ml, Flange diam. 138 mm, Body diam. 140 mm, Height 270 mm, Pressure 1.5* bar, DN 100	1 pz
LLG09142076	Reaction flasks, flat flange, DURAN®, 6000 ml, 7320 ml, Flange diam. 184 mm, Body diam. 215 mm, Height 320 mm, Pressure 1.0* bar, DN 150	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

