

**Lithium batteries, Energizer****Voltronic****Lithium batteries, Energizer®**

- For temperatures from -40 to +60 °C
- Mobile energie on the highest level
- Not rechargeable
- For all high current applications
- 33 % lighter than alkaline AA
- Lasts up to 11 x longer in digital cameras than standard alkaline batteries
- 20 years shelf life



Code	Description	Packaging
LLG09012929	Lithium batteries, Energizer®, Type FR03/L92/AAA/Micro, Volt 1.5	2 pz
LLG09012928	Lithium batteries, Energizer®, Type FR6/L91/AA/Mignon, Volt 1.5	2 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

