



Accessories for Shaker MS 3

IKA

Accessories for Shaker MS 3



Code	Description	Packaging
LLG09729943	Accessories for Shaker MS 3, Type MS 3.4, Microtitre plate attachment	1 pz.
LLG09729945	Accessories for Shaker MS 3, Type MS 1.21, One-hand insert	1 pz.
LLG06237469	Accessories for Shaker MS 3, Type MS 3.5, PCR attachment	1 pz.
LLG09729942	Accessories for Shaker MS 3, Type MS 3.1, Standard attachment	1 pz.
LLG09729949	Accessories for Shaker MS 3, Type MS 1.34, Test tube insert, plain, for user to drill	1 pz.
LLG09729946	Accessories for Shaker MS 3, Type MS 1.31, Test tube insert, with 14 holes for tubes diameter 10mm	1 pz.
LLG09729948	Accessories for Shaker MS 3, Type MS 1.33, Test tube insert, with 4 holes for tubes diameter 16mm	1 pz.
LLG09729947	Accessories for Shaker MS 3, Type MS 1.32, Test tube insert, with 6 holes for tubes diameter 12mm	1 pz.
LLG09729944	Accessories for Shaker MS 3, Type MS 3.3, Universal attachment	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

