



Turbine stirring rotors, 3-blade with guide ring, stainless steel 1.4305

Turbine stirring rotors, 3-blade with guide ring, stainless steel 1.4305

The 3-blade turbine stirring rotor produces axial flow; the guide ring surrounding the blades reduces shear forces, and offers protection when contacting the wall or bottom of the vessel. Recommended for medium to high speed operations.



Code	Description	Packaging
LLG09197030	Turbine stirring rotors, 3-blade with guide ring, stainless steel 1.4305, Ø agitator 30 mm, Ø agitator shaft 8 mm, Length 350 mm	1 pz.
LLG09197031	Turbine stirring rotors, 3-blade with guide ring, stainless steel 1.4305, Ø agitator 50 mm, Ø agitator shaft 8 mm, Length 350 mm	1 pz.
LLG09197032	Turbine stirring rotors, 3-blade with guide ring, stainless steel 1.4305, Ø agitator 70 mm, Ø agitator shaft 10 mm, Length 400 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

